CARNET D'ENTRETIEN
SERVICING & MAINTENANCE MANUAL

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NOTES

CONTENTS

GENERAL REMARKS 24
Intended use 24
Indications and ability 24
Service life 25
Regulation 25
Symbols used in this manual 25
Meaning of symbols appearing on the labels 25

SPECIFICATIONS OF LIGHT ASSIST 3 / LIGHT ASSIST 3 XL 26
General assembly view of the power units 26
Technical specifications 26
Components 27

DRIVING PRACTICES 29
How to drive? 29
How to approach a slope or an obstacle? 29
How to tackle ground sloping sideways? 30
How to brake? 30
How to switch to manual mode? 30

SERVICING AND MAINTENANCE 31
Servicing the equipment 31
Performance checks 31
Waste management 32
Troubleshooting 33

TRANSPORT AND STORAGE 34
Transport 34
Storage 34

SAFETY AND HAZARD INSTRUCTIONS 35
General safety instructions 35
Safety informations about the electrical system 36
Safety recommendations about electromagnetic compatibility 37
Safety information about installation 37
Safety information about driving 37
Safety information about servicing and maintenance 38
Safety information about modifications made to the equipment 38
Important information about the maintenance tools 39
CE marking 39
Important legal notice for the user of this product 39

GUARANTEE 40
Conditions 40
Limits of liability 40

GENERAL REMARKS

Read this servicing & maintenance manual before using this product and keep it for future reference.
The content of this document may change without notice. We recommend you visit our website regularly for any changes or updates.

INTENDED USE

The LIGHT ASSIST 3 range, composed of LIGHT ASSIST and LIGHT ASSIST XL, aid to pushing and braking the manual wheelchair that should only be used by the person accompanying the user. This lightweight power add-on transforms a manual wheelchair into a power wheelchair.

For safety reasons, the LIGHT ASSIST range should only be used by people who :

- assist the user,
- have been told how to handle it,
- are mentally and physically able to use a control unit to operate the motorised wheelchair and its functions.

Users and his attendant must follow the general safety and use operating instructions. As the manufacturer, BENOIT SYSTEMES reserves the right to change product specifications without notice.

Please avoid using the LIGHT ASSIST range before receiving instructions about the equipment from an approved dealer.

For any question about the use, servicing or safety of your device, contact your dealer.
If you have further questions, you can contact:
BENOIT SYSTEMES, 7 rue du Pont, 21450 Billy-lès-Chanceaux, FRANCE
Tel : +33 3 80 96 51 25
export@benoitsystemes.com

INDICATIONS AND ABILITY

The LIGHT ASSIST 3 range is compatible with most manual and tilt-in-space wheelchairs (see the compatibility guide on ourwebsite). It is suitable for caregivers accompanying people who are unable to walk or have limited mobility due to:

- Paralysis,
- Loss of a limb (amputation of the leg),
- A limb problem or deformity,
- Spasms / injuries to the joints,
- Strokes and brain injuries,
- Neurological disorders (multiple sclerosis, Parkinson's, etc.),
- Illness, such as heart and circulatory disease, balance disorder, as well as elderly people with enough strength in the upper body,

Caregivers are individuals who are mentally and physically capable of using a control unit to safely operate the device and its functions.

The prescription conditions must consider all factors for which the LIGHT ASSIST 3 range provides a relevant answer. A PM&R physician (Physical and Rehabilitation Medicine), assisted by an occupational therapist and a physiotherapist, will assess the suitability of the adaptive comfort solution with regard to the patient's condition and abilities, environment, and physical and cognitive capacities.

A validation test will be carried out in the presence of the multidisciplinary team and the patient. The purpose of this test will be to confirm the best match between :

- The physical and cognitive capacities of the patient, which should enable him to safely handle the device.
- And the device itself to meet the need of this particular disability compensation.

SERVICE LIFE

BENOIT SYSTEMES considers a service life of its products at 5 years, as long as the product is used for its intended purpose and obeying use, safety, and maintenance instructions. However, this service life can be exceeded when the product is serviced and maintained carefully and used cautiously, appropriate to its technical specifications
The service life may also be considerably shortened following extreme or inappropriate use. The service life set by BENOIT SYSTEMES is for information only and does not constitute an equivalent or additional contractual guarantee.

REGULATION

The LIGHT ASSIST range, designed by BENOIT SYSTEMES, is considered a "Medical Device" in accordance with Regulation (EU) 2017/745 of the European Parliament and of the Council of 5 April 2017 on medical devices, amending Directive 2001/83/EC, Regulations (EC) No 178/2002 and No 1223/2009, and repealing Council Directives 90/385/EEC and 93/42/EEC.

SYMBOLS USED IN THIS MANUAL

In this manual, hazard statements are indicated by symbols. These symbols are accompanied by a warning indicating the risk level.


WARNING
Indicates a dangerous situation that, if not avoided, is likely to cause severe or even fatal injuries.

CAUTION
Indicates a hazardous situation that, if not avoided, is likely to cause in minor or slight injuries.

INFORMATION
Gives advice, recommendations and useful information for efficient, problem-free use.

MEANING OF SYMBOLS APPEARING ON THE LABELS


	Manufacturer
MD	Medical device
CE	CE marking in compliance with Annex V of European Regulation 2017/745/EEC
SN	Power add-on serial number. It is located under the power add-on, between the battery brackets and on the first page of the servicing & maintenance manual.
	Read the instructions for use carefully - indicates the need for the user instructions and safety guidance carefully.
	Keep dry
	Maximum speed
	Maximum user weight
	This product contains parts that can be recycled and substances potentially harmful to the environment. It must be directed to appropriate waste treatment streams (see "Waste management," page 32).
	UDI (Unique Device Identification) number enabling traceability of medical devices.

SPECIFICATIONS OF LIGHT ASSIST 3 /LIGHT ASSIST 3 XL

WARNING
It is prohibited to modify the equipments without the producer's authorisation.

GENERAL ASSEMBLY VIEW OF THE POWER UNITS

TECHNICAL SPECIFICATIONS

	LIGHT ASSIST 3	LIGHT ASSIST 3 XL
Weight of the power unit (kg)	7,5	9,0
Min. dimensions (mm)	**Width :** adjusts to the width of the wheelchair without extending beyond the handrims. **Depth :** 280 **Height :** 270	
Max. user weight* (kg)	125	150
Max. speed (km/h)	6	6
Max. safe slope (°)	3	3
Max. crossable obstacle (mm)	50	50
Max. stopping distance ** (m) on level and sloped surfaces	1.0 Compliant with EN 12184 standard	1.0 Compliant with EN 12184 standard
Suitable wheelchairs***	Compatible with most manual wheelchairs (Class A wheelchairs)	

** If you exceed this limit, you might damage the frame or fastenings of your equipment, which could cause serious injuries for you or someone else. Ensure the combined weight of the user and the device does not exceed the allowable load of the manual wheelchair (refer to the manual wheelchair's technical specifications).*

*** Stopping distance varies depending on the gradient. It is shorter when driving uphill, but slightly longer when driving downhill.*

**** A compatibility guide is available on our website to know the wheelchair models on which the Light Assist 3 can be installed, in compliance with the general safety and performance requirements set out in European Regulation 2017/745 on medical devices.*

COMPONENTS

The power add-on consists of several components :

1. Power unit
2. Roller covers
3. Rollers
4. Axles
5. Handle
6. Carrying and unlocking handle
7. Control handle
8. Static handle
9. Connector to charge the battery
10. Anti-tip
11. Pin locker
12. Battery
13. Battery handle
14. Power connector (red)

POWER UNIT (1)

It consists of two electric motors equipped with rollers (3), clutch mechanisms, and the battery support.

POWER SUPPLY

The power supply consists of a 282 or 484 Wh lithium-ion battery (12), installed on a bracket that enables automatic installation and connection to the power unit.

It is delivered charged and turned off, except when transported by air.

We decline all liability for damage caused by using other batteries not supplied by us.

A battery user manual is supplied with the product to assist with commissioning and maintenance.


CONTROL HANDLES (7)

All the travel functions of the LIGHT ASSIST are operated by the control handles. Two handles are supplied with the power add-on :


- a throttle handle
- a static handle

The handles are positioned to right or left on the rear handles of the wheelchair or using clamps attached to the stiffener bar or U-shaped push handle bar.

THROTTLE HANDLE

STATIC HANDLE

It is positioned on the opposite side of the throttle handle and completes the control system.

ANTI-TIP (10)

It must always be set up, wheel bottom facing. It has been adjusted by your technician to ensure safe driving while allowing obstacles to be crossed up to 5 cm.

⚠ BENOIT SYSTEMES waives all liability if an accident occurs while the anti-tip is not fitted.

ACCESSOIRES (OPTIONAL)

BENOIT SYSTEMES offers various optional accessories such as:

- **Mudguards, an essential accessory for outdoor trips in bad weather.**
- **A car charger with cigarette-lighter plug,** convenient for charging the battery when traveling.

Only original BENOIT SYSTEMES accessories must be used.

For more information on accessories, please contact your dealer.

ⓘ INFORMATION
The device is programmed at the factory during manufacturing. Any modification to the programming may only be carried out by a BENOIT SYSTEMES technician or a person trained by BENOIT SYSTEMES.

DRIVING PRACTICES


Before making your first journey, you should familiarise yourself with operation of the motorised wheelchair and all the accessories. Take your time to test all the functions and driving modes.

PROPERLY SEATED = SAFE DRIVING

Before each journey, check that:

- The user is positioned safely in their chair (seat belt if needed).
- The battery is sufficiently charged for the journey you're planning to make.
- The safety anti-tip is clipped into the rear handle. Check the mounting direction (small wheel downwards) and that the locking pin is securely fitted.
- Your tires are correctly inflated (recommended pressure 8 bar).

WARNING
Do not use any other motorised device or accessory at the same time as the power add-on. If you wish to use another device, please remove the power add-on from the wheelchair for your own safety.

HOW TO DRIVE?

Press **Start / Stop** [1] button to set LIGHT ASSIST in motion.

Turn progressively the throttle handle to set the wheelchair in motion, after selectiong the desired direction of travel (**Forwards/Backwards button** [2]). The speed is proportional to the travel of the handle. Use the **limiter** [3] if needed, to decrease or increase the device's maximum speed.

For a pleasant journey, adjust the speed to suit your own pace.

When the throttle handle is released, it automatically returns to the neutral position, the wheelchair stops moving, and the parking brakes engage.

HOW TO APPROACH A SLOPE OR OBSTACLE?

The power add-on is designed for safe driving within the slope and obstacle limits stated in the technical specifications.

Approach climbs with sufficient speed and a fully charged battery. Effort can cause a voltage drop leading to progressive stop due to an insufficient charge level.

On a hill start, turn the control handle progressively to avoid abrupt movement of the wheelchair, which can cause rearward imbalance.

When the slope is steady and the speed remains sufficient to keep moving forward, you can continue without adjusting the settings. If the slope changes and the speed becomes too slow, you can help the wheelchair climb by pushing manually, or use the limiter to slightly increase the maximum speed until it stabilizes.

Reduce speed downhill. If the battery is fully charged and speed is excessive, an overvoltage may occur and cause a sudden stop.

Approach any obstacle perpendicularly, preferably in reverse, so both wheels climb simultaneously. First ensure its height is compatible with the anti-tip or the wheelchair's capability.

Approach the obstacle slowly, gently tilt the wheelchair backward to take the weight off the front wheels, and continue moving forward.

If the caster wheels cannot be positioned on the obstacle, or if crossing it is not possible in forward motion, perform the maneuver in reverse, with or without momentum depending on the height to be crossed.

⚠ WARNING
BENOIT SYSTEMES accepts no liability for accidents resulting from incorrect installation or positioning of the anti-tip.

HOW TO TACKLE GROUND SLOPING SIDEWAYS?

When the ground slopes sideways, the wheelchair is drawn naturally towards the lower side. For this reason, you should correct the direction.

The LIGHT ASSIST is steered using arm strength, by applying lateral pressure on the control handles, like a manual wheelchair. Use the handles to guide the wheelchair, anticipating movements as much as possible, especially on slopes.

Reduce speed downhill.

Never traverse steep slopes sideways to avoid the risk of the wheelchair stopping or lateral tipping.

HOW TO BRAKE?

To prevent any uncontrolled movement of the motorised wheelchair, the motors are fitted with automatic brakes that are applied under the following conditions:

- The start/stop switch is in the 'stop' position.
- The start/stop switch is in the 'start' position and the control handle released.
- Power supply is interrupted.

The brake is released automatically as soon as the wheelchair starts moving, when the handle is operated. The safest and simplest way to stop your equipment is to release the throttle handle, in order to control the wheelchair'stopping.

In the event of a power supply fault, the brakes are operated automatically.

If the controller or throttle handle fails without a power failure, the wheelchair first stops (as if the control handle was released), then the brakes are operated automatically

When you need to push the wheelchair for precise maneuvers or in case of emergency, the power add-on can be disengaged.

HOW TO SWITCH TO MANUAL MODE?


If needed, the wheelchair can be used manually without removing the power add-on. Turn the handles to disengage the device

Follow these steps :

Disengaged handles position

1 *Extend the handles.*
2 *Turn the handles outward, then pivot them downard*

Telescopic handles

SERVICING AND MAINTENACE

⚠ CAUTION
It is important to follow the servicing and maintenance programme for the equipment to keep it in good condition.

WARNING
If you're not sure you can service your equipment yourself, contact your trained dealer.

SERVICING THE EQUIPMENT

- Before undertaking any cleaning, ensure the equipment is turned off.
- Clean the rollers and their housing regularly with a screwdriver.
- Once a week, run a vacuum cleaner around the motor, to remove all accumulated dust.
- Clean the entire power unit with a damp (not wet) cloth.
- Clean the electrical components and control unit with clean and dry cloth.
- Never use abrasive or alcohol-based detergents.
- Wipe the device immediately if wet or damp after cleaning or use.

WARNING
Never wash your wheelchair using a hose or pressure washer.

Only use common household cleaning and disinfection products, except abrasive or alcohol-based detergents.

WARNING
Recommended tyre pressure is 8 bar. It must never be below 6 bar.

Check tyre pressure regularly to ensure that it is always as close as possible to the recommended value. Use only the tyre brands and models recommended by BENOIT SYSTEMES and fitted at installation (SCHWALBE Marathon Plus or PRIMO Sentinel). Contact your dealer if you cannot find them.

Never replace the pneumatic tyres with solid tyres.

PERFORMANCE CHECKS

After having cleaned or repaired the device, you must always check it is working correctly before using it:

- Turn on the throttle handle.
- Test the braking system on a flat surface with clear space of at least 1 m around the wheelchair : turn on the throttle handle and select the desired direction. Turn slowly the handle and move forward or backward until you hear the parking brakes are disengaged (it is possible that the wheelchair will move forward). Release immediately the handle. The brakes must then operate and stop the wheelchair.
- Repeat the test a further three time, forward and backward.
- Activate all the electrical options to confirm they are working correctly.
- Use the device in all driving modes to confirm that it operates as before.
- Check the condition of the clutch handle ball joints and ensure they move freely without binding.

⚠ CAUTION
Regularly check device condition to ensure safe use. If in doubt about the performance of your device, contact your dealer.

It is recommended that the items shown in the table are inspected at the stated intervals.
Your approved dealer should repair or replace any loose, worn, distorted or damaged parts before using your device

Servicing and inspection programme	Daily*	Monthly	Quaterly	Annually
Check the battery charge indicator and the charger	●			
Check that the control handles and its components are neither distorted nor damaged	●			
Check that all the parts and devices are firmly attached	●			
Test the braking system	●			
Inspect the condition of the tyres and wheel fastenings	●			
Check that the cables and connectors are in good condition, are not tangled and are not in the way	●			
Inspect the condition of the seat you use		●		
Inspect all the bolts, nuts and fastenings to confirm they are not damaged or loose		●		
Check that the battery terminals are tightened			●	
Check the condition of the clutch handle ball joints and ensure they move freely			●	
Perform a full inspection, safety check and service by a BENOIT SYSTEMES approved dealer				●

* Before each use

Correct servicing will improve your vehicle's performance and extend its service life, as well as reducing the risk of an accident.

WASTE MANAGEMENT

Disposal and recycling must comply with applicable national and local regulations.

Protect the environment by recycling this product via appropriate channels:

- Power add-on and accessories (Waste Electrical and Electronic Equipment): to be disposed of at a waste collection center or to your dealer.
- Battery: to be disposed at waste collections centre.
- Packaging and packing materials: to be disposed of at a waste collection centre.

TROUBLESHOOTING

The power add-on doesn't work	
Possible causes	**Solutions**
Battery incorrectly locked in place	Remove the battery slightly and re-engage it
Throttle handle wire poorly connected	Check the connection
Red connectors damaged or loose	Check that the connectors are in good condition and that the contacts are correctly positioned
Battery fault	Contact your dealer
A connection is faulty	Check terminals and connections
The throttle handle switches on normally but the power add-on doesn't work	
Possible causes	**Solutions**
The power add-on is in the disengaged position	Engage the two handles
The power setting is too low	Increase the power
Rapid wear of the tyre treads	
Possible causes	**Solutions**
Poor quality tyres	Fit recommended tyres
Tyres under-flated	Must not be below 6 bar
Insufficient pressure of the rollers on the tyres	Contact your dealer
Tyre worn on the inner wall	
Possible causes	**Solutions**
The roller cover is touching the tyre	Contact your dealer
The bracket for the wheelchair wheel is loose	Tighten this bracket correctly
The wheelchair goes in a circle	
Possible causes	**Solutions**
A manual brake of the wheelchair is on	Release it
A motor is faulty	Contact your dealer
The power add-on is making noise	
Possible causes	**Solutions**
Part of the wheelchair is touching a wheel	Visually check the condition of the wheelchair
A motor is faulty	Contact your dealer
A bearing is faulty	Contact your dealer
A roller cover is touching the tyre	Contact your dealer
A clutch handle is not staying engaged	
Possible causes	**Solutions**
Incorrect adjustment of the ball joints	Contact your dealer
Distortion (following an impact for example)	Contact your dealer
The power add-on does not unlock properly	
Possible causes	**Solutions**
Clutch wire misadjusted	Contact your dealer

This diagnosis is only for guidance. Once the fault code is seen, please contact your BENOIT SYSTEMES trained dealer.

TRANSPORT AND STORAGE

TRANSPORT

By air

⚠ WARNING
Before departure, contact the airline for all specific conditions and instructions to be followed.
Check that removable parts are firmly secured or packed separately and labeled to avoid loss during loading/unloading. Take this manual with you.

By car

For safety reasons, the power add-on will be put in the boot of the car. The device must be removed when the manual wheelchair is used as a seat in a car. BENOIT SYSTEMES offers an optional charger for cigarette-lighter sockets.

⚠ WARNING
BENOIT SYSTEMES and it representatives refuses to accept any liability in case of accidents and/or damage on the wheelchair resulting from inappropriate use or failure to follow these instructions and/or applicable national regulations.

STORAGE

If the device is not used for more than a week, you should prepare it for storage in order to extend its service life and the battery lifetime.

Follow these instructions :

- Store the device at a temperature of 15°C and avoid extreme temperatures in order to extend the service life of the device and batteries.
- The components are tested and approved for the temperature ranges below :
 * The range of temperatures permitted for storing the device is between **-40° and 65°C.**
 * The range of temperatures permitted for storing the batteries is between **-20° and 50°C.**

WARNING
Never store your equipment

- Outdoors.
- In direct light.
- Close to a source of direct hear.
- In a damp environment.
- In a cold environment.
- With the battery connected (even if the control handle is turned off).

⚠ CAUTION
After an extended storage period at low temperature (–25°C), it is possible for the power add-on to turn itself off after being turned on; this situation may persist for about 2 minutes. In this case, place it at room temperature, wait for a few minutes, then try turning it on again

For battery transport and storage conditions, refer to the battery user manual supplied.

SAFETY AND HAZARD INSTRUCTIONS

GENERAL SAFETY INSTRUCTIONS

WARNING
Risk of damage to the equipment or serious, even fatal, injury if the control handles are used incorrectly.

- Undertake a test session with a healthcare professional before using the control handles.
- If the attendant is not used to operate with control handles, it is recommended to drive slowly and familiarise with the driving and directional characteristics of the wheelchair on a flat, unobstructed route.

WARNING
Risk of damage to the equipment or serious, even fatal, injury.

- Dropping lighted cigarettes onto a padded seating system can start a fire likely to cause damage to the equipment or serious, even fatal, injuries.
- The occupants of a motorised wheelchair are more particularly exposed to a risk of serious or even fatal injury associated with such fires and the resulting smoke, as they are not always able to get away from the wheelchair – Do not smoke when using this device.

WARNING
Risk of damage to the equipment or serious, even fatal, injury.

- Inappropriate supervision or servicing are likely to result in damage to the equipment or serious, even fatal, injuries, following breathlessness or ingestion of parts or materials.
- It is essential for children, pets or people with physical/learning disabilities to be supervised closely.
- Don't allow children or anyone else to use your device.

WARNING
Risk of injury if the equipment is used for a purpose other than that described in this manual.

- The device must be used in accordance with the manufacturer's instructions. It must not be used for a purpose other than that for which it is designed. Pay special attention to the safety instructions.

WARNING
Risk of damage or injury if the device is started unintentionally.

- Do not sit on the motorised wheelchair during the adjustment phase.
- Disconnect the power supply from the device before getting onto it, getting off it or handling objects.
- Do not push the wheelchair by hand when the brakes are applied and when the power supply is disconnected.

WARNING
Risk of falling from the motorised wheelchair.

- Do not move forward on the seat, do not lean forward and do not stretch backwards beyond the top of the seat-back, such as to reach an object.

WARNING
Risk of injury if driving under the influence of drugs or alcohol.

- Never drive if you are under the influence of drugs or alcohol. If necessary, ask for help from someone else with appropriate physical and mental abilities.

WARNING
Risk of injury if the contact is broken while the motorised wheelchair is moving, leading to the wheelchair stopping suddenly and violently.

- Never operate the stop button nor disconnect a cable while the wheelchair is moving, except in an emergency.
- If you have to brake in an emergency, simply release the control handle to stop the wheelchair. Apply the wheelchair brakes on.

WARNING
Risk of damage to the equipment or serious injury due to heat exposure.

- Storing or using the device close to a naked flame or flammable products can result in damage to the equipment or serious injuries.
- To ensure the equipment is reliable, avoid exposing it to extreme conditions, heat or humidity.
- Never store your device outdoors.

⚠ **CAUTION**
Risk of injury if the maximum permitted load is exceeded.

- Do not exceed the maximum permitted load (refer to the Technical Specifications section page 6).
- Check that the combined weight of the device and the user does not exceed the maximum permissible load of the manual wheelchair (refer to the technical specifications of the manual wheelchair).
- Never use the device to transport several people.

⚠ **CAUTION**
Risk of injury in the event of inappropriate handling or fall of heavy part.

- During servicing and maintenance operations or when you are lifting certain parts of the equipment, consider the weight of each component, particularly the batteries. Ensure you always adopt the correct posture and don't hesitate to ask for help.

⚠ **CAUTION**
Risk of injury caused by moving parts.

- Ensure that the moving parts of the motorised wheelchair, such as the wheels or control handles, do not cause injuries, particularly when children are present.

⚠ **CAUTION**
Risk of fire or fault when connecting electrical equipment.

- Do not connect your motorised wheelchair to any equipment that has not been explicitly permitted by BENOIT SYSTEMES. Rely on your trained dealer for all electrical installations.
- Do not use any other motorised device or accessory at the same time as power add-on. If you wish to use another device, please remove the power add-on from the wheelchair for your own safety.

SAFETY INFORMATION ABOUT THE ELECTRICAL SYSTEM

WARNING
Risk of damage to the equipment or serious, even fatal, injury due to fire.
Inappropriate use of the motorised wheelchair can result in the emission of smoke or sparks, and heat may be given off.

- The device must not be used for a purpose other than that for which it is designed.
- If the device begins to emit smoke, sparks or heat, stop using it and have it checked immediately.

WARNING
Risk of injury or damage associated with shortcircuits.
The pins of electrical connectors may remain live even when the system is switched off.

- Check that the cable insulation is sound and that the connectors are correctly plugged in, so that no live parts are accessible or in contact with conductive materials.
- Insulate the pins of the connectors when they are disconnected.

WARNING
Risk of damage to the equipment or serious, even fatal injury due to moisture exposure or liquids, or corrosion of electrical components.

- Reduce the exposure of electrical components to water and/or liquids.
- Electrical components damaged by corrosion must be replaced immediately.
- The power add-ons frequently exposed to water/ liquids can require electrical components to be replaced more frequently.

WARNING
Risk of serious, even fatal injury due to electric exposure.
Failure to take note of these warnings can cause an electrical discharge resulting in serious or fatal injuries, or damage to the electrical system.

- The (+) battery terminal (red or brown cable depending on the model) must be connected to the red (+) power supply cable on the control box.
- The (-) battery terminal (black or blue cable depending on the model) must be connected to the black (-) power supply cable on the electronics unit.
- Ensure to never connect the two poles of the battery (direct contact, tools or conductive material).
- Replace immediately cables with damaged insulation, frayed cables or cut cords.
- Do not disassemble or remove any battery connectors, covers or fuses.

⚠ **CAUTION**
Risk of damage to the equipment due to electrical system failure.
A fault in the electrical system can cause unusual behaviour.

- In the event of a fault, turn the control unit off and back on again.
- If the problem persists, disconnect or remove the battery.
- If in doubt or for any question, don't hesitate to contact your dealer.

SAFETY RECOMMENDATIONS ABOUT ELECTROMAGNETIC COMPATIBILITY.

Electromagnetic fields, may have an influence on the operation of the vehicle. The electronic equipment used in this device can also cause weak electromagnetic disturbances, however these are below the legal limit. Therefore please consider the following comments:

WARNING
Risk of incorrect operation following electromagnetic emissions.

- Do not use a portable transmitter nor communication devices while the device is in use (radio devices or mobile phones).
- Avoid going near powerful radio or television transmitters.
- If your wheelchair should start moving by itself, or if the brakes happen to be released, turn the device off.
- Report to your dealer all undesirable movements of the wheelchair that occur, even the electric brakes being released.
- Do not add any electrical or electronic accessories to the device.

SAFETY INFORMATION ABOUT INSTALLATION

WARNING
Risk of damage to the equipment or serious, even fatal, injury if improperly adjusted.
Prolonged use of the device with settings not complying with the specifications risk causing the vehicle to become unstable and causing damage to the equipment or serious, even fatal, injuries.

- Performance adjustments must only be made by people trained by BENOIT SYSTEMES who fully understand this process and the user's capabilities.
- After adjusting the device, ensure that it operates according to the parameters defined during the installation procedure.

If the device does not operate correctly, contact your dealer.

WARNING
Risk of damage to the equipment or serious, even fatal, injury de to faulty mounting.
Inadequately tightened or missing fastenings can cause the vehicle to become unstable, resulting in damage to the equipment or serious, even fatal, injuries.

- After any adjustment, repair or servicing and before use, ensure that all the fastenings are fitted and suitably tightened.

⚠ **CAUTION**
Risk of damage to the equipment or serious injury due to faulty handles.

- Ensure that the control unit is securely mounted and that all screws are correctly tightened after adjustment.
- In the event of faulty control unit, immediately turn off the device's electronic system.

⚠ **WARNING**
Pinch points can arise when making adjustments.
Be careful of your fingers.

⚠ **CAUTION**
Check all cables.
Damage caused to the cables of the wheelchair can cause safety problems and must be reported immediately to a trained technician.

SAFETY INFORMATION ABOUT DRIVING

WARNING
Risque of damage to the equipment or serious, even fatal, injury.

- Adapt your speed to the environment in which you are travelling. To control handling of the equipment, it is better to start on flat ground.
- Always ensure you can be seen from afar, especially if you use your motorised wheelchair at night.
- Pay close attention to traffic on the road.
- Don't forget that car or lorry drivers do not expect to see a wheelchair running at a speed greater than that achieved by propulsion. If there is the least danger, don't take the risk of crossing the road.

- The device can stop operating suddenly if a defect operates : we recommend using the seat belt when using the power add-on.
- Ensure to adjust the speed appropriately to safely overcome obstacles.
- Any malfunction of the control unit is liable to result in unexpected/erratic movements leading to damage to the equipment or seriouss, even fatal, injuries.
- If an unexpected/erratic movement occurs, stop using immediately the device and contact your dealer.
- Control handles must be replaced if it is damaged or weakened.

WARNING
Risk of injury if the wheelchair tips over.
- Refer to the technical specifications (maximum safe slope).
- As far as possible, avoid driving on slippery, icy or oily surfaces when there is a risk of losing control of the vehicle. If you cannot avoid travelling over such surfaces, drive slowly and with the greatest care.
- Ensure that there are no objects in your path likely to damage the mechanism of your device.
- Do not cross an obstacle while moving forward.
- Avoid passing over thresholds or obstacles above the specified maximum height to not to threaten your stability. Remember to arrange for a slope, if necessary.
- Never try to go down a series of steps with your device.
- Avoid making sudden movements of the control unit od sudden changes of direction while the motorised vehicle is moving.
- Never use the motorised wheelchair to transport several people.
- Check that the user weight combined with the device do not exceed the maximum total permitted load (refer to the technical specifications of the wheelchair and the power add-on).
- Remember that the device brakes or accelerates when you change the speed of travel while the vehicle is moving.
- Never remove the anti-tip before or during the use of the device.
- **If your wheelchair is a standing model (e.g., LEVO Summit or equivalent)**: never use the power add-on when the wheelchair is in the standing position. Always remove the control unit before moving into the upright position.

WARNING
Risk of injury if you collide with an obstacle while passing through narrow passageways.
Pass through narrow passageways at the slowest speed of travel and with the greatest care.

SAFETY INFORMATION ABOUT SERVICING AND MAINTENANCE

WARNING
Risk of damage to the equipment or serious, even fatal, injury.
Incorrect repair or maintenance of this device by users/carers or by unqualified technicians are liable to result in damage to the equipment or serious, even fatal, injuries.
- Do not undertake any maintenance task other than those described in this servicing manual. These repairs and/or servicing must be carried out by a qualified technician. Contact a BENOIT SYSTMES dealer or technician.

⚠ CAUTION
Risk of accident if incorrect maintenance performed.
- For safety reasons and to avoid accidents related to passive wear, it is vital that this device should be checked once a year under normal use conditions.

SAFETY INFORMATION ABOUT MODIFICATIONS MADE TO THE EQUIPMENT

⚠ CAUTION
Using unapproved components or accessories is liable to result in injuries or to damage the motorised vehicle. Using accessories and parts not approved by BENOIT SYSTEMES may alter the stability of the motorised wheelchair and increase the risk of tipping over.
- Use only accessories and parts approved by BENOIT SYSTEMES for this equipment.
- Electrical and electronic components not approved by BENOIT SYSTEMES for this device are liable to cause fires and electromagnetic damage.
- Use only electrical and electronic components approved by BENOIT SYSTEMES for this equipment.
- Ensure you always quote the serial number of the equipment when you order spare parts.

IMPORTANT INFORMATION ABOUT THE MAINTENANCE TOOLS

Certain maintenance tasks described in the servicing manual can be performed by the user without problems require the use of appropriate tools. If you don't have these tools, we advise you not to perform these maintenance tasks. If necessary, we recommend contacting your dealer.

CE MARKING

Conformity was assessed and CE marking applied in compliance with European Regulation 2017/745/EEC and only applies to the complete product.
CE marking is invalidated if components or accessories are replaced or added without having been approved by BENOIT SYSTEMES for this product.

IMPORTANT LEGAL NOTICE FOR THE USER OF THIS PRODUCT

Incidents due to device malfunctions and resulting in serious physical injuries must be reported to the producer and to the competent authority of the Member State where the user is located;

IN FRANCE
- Manifacturer : BENOIT SYSTEMES, 7 rue du Pont 21450 BILLY-LÈS-CHANCEAUX
- Competent authority : ANSM, 143 Boulevard Anatole France 93200 SAINT DENIS

IN THE MEMBER STATES OF THE EUROPEAN UNION
- Your distributor.
- National authorities responsible for incidents with medical devices.

You will find an overview of the responsible authorities by following the link : http://ec.europa.eu/growth/sectors/medical-devices/contacts/

GUARANTEE

The guarantee period for the power add-on, battery and charger is 24 months from the invoice date.

All spare parts come with a 6-month warranty from the date of delivery, while all batteries purchased separately are covered by a one-year warranty.

To claim against the guarantee, just make a note of the serial number of your equipment shown on the first page of this manual and contact your approved dealer or to register on our website. The serial number is also under the power add-on, between the battery brackets.

CONDITIONS

The guarantee is only valid if the product is used under the specified conditions and for the use for which it is designed.

The guarantee may not be applied if the product requires repair or replacement for one of the following reasons:

- The power add-on, or one of its parts, has not been maintained according to the manufacturer's recommendations, as stated in the servicing manual or parts other than those specified for this equipment are used.
- The power add-on, or one of its parts, was damaged as a result of negligence, an accident or inappropriate use.
- The power add-on, or one of its parts, has been converted or a repair made before having notified the after-sales service department.

This guarantee is limited to the free-of-charge repair or replacement of parts recognised as faulty in the manufacturer's workshops or those of its representative (including labour and carriage). This repair may be performed using new parts or standard exchange parts, depending on the assessment made by the manufacturer or its representative.

Work carried out under the guarantee does not extend the guarantee.

LIMITS OF LIABILITY

BENOIT SYSTEMES declines all liability for damaged related to:

- Failure to comply with the user manual.
- Incorrect or improper use of the equipment.
- Normal wear and tear.
- Inappropriate storage of the equipment.
- Incorrect assembly or dismantling by the user or other parties.
- Technical modifications on and in the equipment.
- Unauthorised modifications and/or use of spare parts not manufactured or recommended by BENOIT SYSTEMES.

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This medical device is a regulated health product and bears the CE marking. As manufacturer, BENOIT SYSTEMES declares that this products is Class I medical devices and complies with EU Regulation 2017/745/CEE on medical device.

7 Rue du Pont, 21450 Billy-lès-Chanceaux FRANCE | (+33) 3 80 96 51 25 | export@benoitsystemes.com

WWW.BENOITSYSTEMES.COM