

**WELL-BEING ACCESSORIES** 



Follow us!









## KITCOOL 2.0


Kitcool 2.0 is a detachable and ergonomic accessory for manual wheelchairs that enables the user to adopt the perfect resting position whenever they choose, by simply tilting backwards without having to transfer into a different seat.

Foldable, lightweight and compact, you can take it everywhere and store it easily.





# HOW DO YOU SET IT UP?





LET'S RELAX!

# ADVANTAGES OF KITCOOL 2.0



# **F**(1)



#### HEALTH

Prevents and helps heal pressure sores, improve circulation, can alleviate postural fatigue, eases breathing and swallowing.

#### WELL-BEING

Helps the user rest and relax as they can change position as easily as an able-bodied person changes chair. INSTANT
INSTALLATION
WITHOUT TOOLS



#### SECURITY

Extremely stable, reliable and robust.



#### EASY SET UP

of the tilt and the headrest.



#### FOLDING, LIGHTWEIGHT AND SLIM LINED

It travels smoothly by car and plane, and can be stored easily.



#### TRAINING AID

Use it to learn to do 'wheelies' with no risk of falling backwards.



Find a relaxing position by yourself or let yourself be helped.

Provides mutual comfort for the carer and the user.

Safe assistance to perform a wheelstand.



## TECHNICAL **DETAILS**

Size of folded Kitcool 2.0 — L600mm x W460cm x H100mm Cushion Veltosuede - 100
Tilt angle from 25° to 60° Veltosuede - 100% polyester

Weight — - 3.7 kg Warranty — 2 years

#### IS IT COMPATIBLE WITH MY WHEELCHAIR?

Kitcool 2.0 can fit on all manual wheelchairs with rear handles.

Contact us!

#### 2 MODELS **OF KITCOOL 2.0**

- For seat width from 26 cm to 35 cm
- For seat width from 36 cm to 52 cm

min 9 cm Required handles length -Required handles diameter - from 20 to 35 mm above the backrest Push handles position from 26 cm to 52 cm Seat width -



## NEW IN JUNE 2022



No tools required for installation



Simple adjustment of the headrest cushion



Simplified mounting system



Stored headrest adjustment





# HOW DO YOU SET IT UP?

# YOUR WHEELCHAIR SET UP THE BRACKET ON THE WHEELCHAIR SET UP THE BRACKET ON THE WHEELCHAIR SET UP THE BRACKET ON THE STOP SLIDE THE METAL BAR OF THE TRAY INTO THE BRACKET THE WHEELCHAIR

# ADVANTAGES OF ACTI'TAB









REMOVABLE AND EASILY TRANSPORTABLE ADAPTS TO ALL MANUAL AND POWER WHEELCHAIRS (contact us) EASY INSTALLATION, NO TOOLS REQUIRED FULLY ADJUSTABLE (HEIGHT AND TILT)



FLIP DOWN POSITION AND ALWAYS AVAILABLE



DISCREET AND SIMPLE



ADJUSTABLE IN HEIGHT AND TILT — LET'S ENJOY!

**Compatible with Light Drive 2.1** 

You can still sit at the table, even if Acti'tab is mounted to your wheelchair.



## TECHNICAL **DETAILS**

3 MODELS **OF ACTI'TAB** 



Tray width 18 cm Tray material -Plexiglas 2 years Warranty



### SEVERAL POSITIONS

## **AVAILABLE**



Flat position for carrying your goods



**Tilted position** for easy reading



Flip down position for easy transfers



#### **ACCESSORIES**

Many accessories are provided including:

- a removable tray with a multi-edge face / a single edge face
- a cup holder
- a non-slip silicone mat...

More information is available on our website.





C E Kitcool 2.0 and Acti'tab comply with the European Regulation 2017/745 on Medical Devices (Class I).